

ergotron[®]

moving you forward

- Helps reduce cart downtime with real-time insights into battery life and more through an online dashboard
- Software communicates directly with the cloud without access to patient records to keep information secure and support HIPAA compliance
- Provides real-time notifications to your IT staff when a battery requires service or attention to eliminate guesswork and unnecessary checks
- Enables IT staff to diagnose issues from a distance to help reduce support time instead of traveling across expansive campuses to locate a cart
- Based on your service level agreement, automated service notifications allow Ergotron to address an event as it occurs to remotely diagnose the situation and dispatch service and parts
- Free to all LiFeKinnex users
- Easy to install via www.ergotron.com/ekinnex without a software download on the workstation; quick to implement with installation per device in as little as one minute

eKinnex[™]
SOFTWARE



FOR
LiFeKinnex[™]

CLOUD-BASED FLEET MANAGEMENT

21st century, patient-centered healthcare needs 21st century technology and equipment. Rethink how you manage your medical carts with a LiFeKinnex Power System with eKinnex[™] technology. This secure, web-based software hosted by Ergotron monitors battery life in real time, helping you proactively identify service and maintenance needs to maintain uptime, freeing up time for your IT team. The HIPPA-compliant portal keeps patient information safe while providing key insights to manage your fleet to meet the needs of fast-paced, modern care.

Detailed power system metrics

Historical power system data, including charge, load and current, lets you quickly identify usage patterns and potential issues.



Dashboard Ergotron

Replace Battery × 1 | Low Battery × 6 | Lost Contact × 5

Carts Ergotron

Filters

Search: Name or Group

Connection: Online Offline

Power Source: Plugged In Battery None

Battery Condition: Good Replace

Charge: 0% — 100% < 100%

Carts (8)

Name	Last Seen	Group	Type	Power Source	Charge	Battery Condition
ERG-TANGENT-PC	2017/01/27 3:30:18 pm	HIMSS17	Unknown	Unknown	0%	Checking
ERGLOANER08	2017/01/24 4:20:51 pm	demo	Unknown	Valence 12RTL 40AH	0%	Checking
NHAZZARD1	2017/02/08 12:06:40 pm	demo	Unknown	Unknown	0%	Checking
SOLSON-TEST	2017/02/08 12:05:57 pm	demo	Unknown	Unknown	81%	Checking
SSQUIRES1	2017/02/02 12:35:44 pm	Unknown	Unknown	Unknown	0%	Checking
SURFACEPC03	2017/02/08 8:51:57 am	demo	Unknown	Unknown	100%	Replace
TTOWNES	2017/02/02 4:36:58 pm	Unknown	Unknown	Unknown	0%	Checking
USER-PC	2017/01/27 3:30:31 pm	HIMSS17	Unknown	Unknown	0%	Checking

Key features:

- Add and remove eKinnex users
- View the state of your hot swap batteries
- Organize carts into groups for easy management
- Sort and filter to identify maintenance needs

Learn more at www.ergotron.com/ekinnex.

At-a-glance insights into your fleet:

- Real-time notifications when a mobile workstation requires service or attention
- Battery status
- Cart hardware information, like last reported location

Download additional resources at ergotron.com.

For more information:

Americas: 800.888.8458 / +1.651.681.7600 / insidesales@ergotron.com

EMEA: +31.33.45.45.600 / info.eu@ergotron.com

APAC: apaccustomerservice@ergotron.com

Custom: custom@ergotron.com

© Ergotron, Inc. 08.05.2023

Literature made in the U.S.A. Content subject to change.

Ergotron devices are not intended to cure, treat, mitigate or prevent any disease.

Patent information available at www.ergotron.com/patents.

WorkFit, StyleView, LearnFit, Neo-Flex, PowerShuttle, LiFeKinnex, CareFit

and eKinnex are trademarks of Ergotron.

