

ergotron®

moving you forward

StyleView 10 Carts

Ergonomic

- 38 cm of effortless worksurface height adjustment lets users sit or stand throughout the day for a healthier way to work

Options

- Three models available to support most laptops and tablets, and all Microsoft Surface devices (see details on back page)

Designed for healthcare

- Small, lightweight and extremely agile
- Sturdy, safe metal construction; powder-coat paint finish on smooth steel surfaces
- Super easy to keep clean
- Includes dual-wheel, 7.6 cm rubber casters; two casters can be locked if necessary
- Includes scanner holder; integrated T-slot accepts additional StyleView accessories and add-on power solutions
- 7-year warranty



LAPTOP
SV10-1100-0



TABLET
SV10-1400-0



SURFACE
SV10-1800-0

WHERE PRODUCTIVITY MEETS STYLE: POINT-OF-NEED CARTS FOR LAPTOP, TABLET AND MICROSOFT® SURFACE™ DEVICES

STYLEVIEW™

StyleView 10 (SV10) carts let caregivers easily move their mobile devices to the point of need. Super petite, lightweight and agile, SV10 carts glide effortlessly across various flooring types and fit nicely into small work areas, like exam rooms and offices.

Designed with ergonomics in mind, carts offer 38 cm of effortless vertical height adjustment so users can choose to sit or stand as they work. Keyboards and displays can be properly positioned for comfortable data entry and screen viewing.

Accessories can easily be added to further enhance functionality and productivity.

STYLEVIEW™



Laptop Cart



Tablet Cart



S-Tablet Cart for Microsoft® Surface™

Part # (Color)	SV10-1100-0	SV10-1400-0	SV10-1800-0
Typical LCD Size*	≤ 17"	≤ 12"	≤ 12"
Worksurface Dimensions	45.7 x 31.8 cm	45.7 x 31.8 cm	58 x 43 cm
Base Footprint	48.3 x 48.3 cm		
Weight Capacity*	≤ 4.5 kg	≤ 1.1 kg	≤ 3.6 kg
Lift	38.1 cm		
Tilt	—	128°	—
VESA	—	MIS-D	—
Cart Weight	15 kg		
Shipping Dimensions	102 x 60 x 60 cm		
Shipping Weight	23 kg		
Warranty	7 years		

* See website for details

SV10 Cart Differentiators

Laptop SV10-1100-0

A bracket holds the laptop in place; inconspicuous design allows easy tilting of the screen without any occlusion of keys or views

Tablet SV10-1400-0

Install tablet in landscape or portrait mode
Tilt the tablet to a preferred input position— work with screen angled vertically, horizontally or anywhere in between. Worksurface provides a convenient platform for a keyboard

Surface SV10-1800-0

Special design of worksurface uses embedded magnets for click-in stability of Surface tablet
Place the Surface in three possible positions: left, right or center
Spacious worksurface holds any Surface tablet with extra area for a mouse or other task-related items



With 38 cm of height adjustability, our patented lift engine technology lets users easily adjust their platform and device to the right ergonomic height, whether they are sitting or standing.

The option to sit and stand throughout the day offers a much healthier way to work!



All SV10 carts include an integrated scanner holder



Integrated hand-brake secures platform height to provide an ultra-stable worksurface



T-slot channel interface on the column allows you to add accessories to your SV10 Cart, such as a Sani-Wipes Holder (97-828) or Storage Bin (97-926-064)

