

# ergotron®

*moving you forward*

## **CareFit Pro Medical Cart**

- Efficient workflows: Customizable drawers, lighting and more streamline daily tasks
- Ergonomic working: Offers full sit-to-stand adjustment for comfortable documentation
- Promotes uptime: Oversee drawer access, track battery life and simplify maintenance with remote fleet management
- Enhanced security: Supports Imprivata identity and access management with optional locking drawers to help keep medications secure
- Professional grade: Designed to last with industry-leading quality testing
- Long-term investment: Adapts to future workflows by adding drawers and accessories



FUTURE-PROOF WORKFLOWS AND PRIORITIZE  
CAREGIVER WELL-BEING


## CAREFIT™

The CareFit Pro Medical Cart sets a new standard in professional-grade design to create home-like environments for clinicians and patients alike. A flexible open architecture system reliably adapts to future and current workflows with ease. Integrated, cloud-based software simplifies fleet management, helps improve uptime and alleviates power concerns.

# CAREFIT™

FULL-FEATURED • HIGHLY CUSTOMIZABLE • INDUSTRY-LEADING

- 1 180° pan for monitor and keyboard
- 2 Touchscreen user interface
- 3 High-capacity configurable drawer system can easily be added
- 4 Infinite brake mechanism enhances patented Constant Force™ Technology for improved usability; electric lift also available
- 5 Integrated 320Wh LiFe battery
- 6 Durable, 12,7 cm dual-wheel casters
- 7 Optional LiFeKinnex™ Power System (245Wh)
- 8 Cantilever monitor mount for focal depth adjustability and full access to worksurface
- 9 LED lighting located on the base, keyboard and worksurface for working in dark environments
- 10 Neutral-angle, dual-grip handle design for improved push/pull/lift ergonomics
- 11 Integrated 7-port USB hub for data and power
- 12 Technology enclosure
- 13 Enclosed cable management
- 14 Flex-platform accessory mounting system
- 15 Promote uptime with eKinnex™ software for proactive fleet management



CareFit Pro with multiple drawer configurations!



	Electric Lift Technology	Constant Force Lift Technology
Part Number	C52-2201-2 (EU) C52-2201-3 (UK, HK, SG) C52-2201-C (CHE)	C52-1201-2 (EU) C52-1201-3 (UK, HK, SG) C52-1201-C (CHE)
Color	White/warm grey	
Typical LCD size	Up to 27" or 9,1 kg	
Worksurface dimensions	35,6 x 40,6 cm, expandable to 35,6 x 58,4 cm	
Base footprint dimensions	50,8 x 53,3 cm	
Lift	68,6 cm: 56 cm + 12,7 cm independent monitor lift	
Warranty	See website	

# CAREFIT™

## CUSTOMIZABLE FOR EVERY WORKFLOW

Need drawer storage?  
If so, how many?

- **One row:**  
Select our cart with patented Constant Force Technology.
- **More than one row:**  
Select our cart with an electric lift.  
Mix and match up to a total of 4 rows (1 Drawer, 1x1) on an electric cart.




1 Drawer (1x1)



2 Drawers (2x1)



5 Drawers (4x1+1)



8 Drawers (4x2)




1 Tall Drawer (1x1)

Supports Imprivata identity and access management.




Save assembly time and reduce your environmental footprint with preconfigured carts.



# CAREFIT™

Configurable drawer system integrated with cart or available as accessories. Configuration examples:



## DO YOU REQUIRE HOT SWAP BATTERY POWER?

Supplement the integrated battery with reliable hot swap power for nearly infinite run times.



## CHOOSE YOUR ACCESSORIES (next page).

Accessories easily snap in place on the left or right of the worksurface without having to remove existing accessories. The direct install T-nut allows you to place your own accessories on the cart.



### Medical Cart with Drawers

Part Number / Availability

1 Drawer (1x1)

Constant Force Lift Technology C52-1211-2 / EU

Constant Force Lift Technology C52-1211-3 / UK

Constant Force Lift Technology C52-1211-C / CHE

1 Drawer (1x1)

Electric Lift Technology C52-2211-2 / EU

Electric Lift Technology C52-2211-3 / UK

Electric Lift Technology C52-2211-C / CHE

2 Drawers (2x1)

Electric Lift Technology C52-22A1-2 / EU

Electric Lift Technology C52-22A1-3 / UK

Electric Lift Technology C52-22A1-C / CHE

5 Drawers (4x1+1)

Electric Lift Technology C52-2251-2 / EU

Electric Lift Technology C52-2251-3 / UK

Electric Lift Technology C52-2251-C / CHE

8 Drawers (4x2)

Electric Lift Technology C52-2281-2 / EU

Electric Lift Technology C52-2281-3 / UK

Electric Lift Technology C52-2281-C / CHE

1 Tall Drawer (1x1)

Electric Lift Technology C52-22B1-2 / EU

Electric Lift Technology C52-22B1-3 / UK

Electric Lift Technology C52-22B1-C / CHE

Accessories

Quad Drawer 98-485

Double Drawer 98-486

Single Drawer 98-487

Double Tall Drawer 98-488

Single Tall Drawer 98-489

# CAREFIT™



## ACCESSORIES



### Document Holder

Replacement for standard cart feature 98-433

Use the kit to quickly attach new accessories to the t-slot without removing other items.



### Direct Install T-Nut Kit

98-432



### Rear Document Holder

Replacement for standard cart feature 98-434



### Bin

98-438

REAR-MOUNT

ACCESSORIES



ACCESSORIES


WORKSURFACE

Supports one worksurface accessory on left and right side



### Slotted Shelf

98-436



### Worksurface Expansion Kit

98-437



### Wipes Holder

98-435



### Supply Tray

98-439

Download additional literature at [ergotron.com](http://ergotron.com)  
 For more information:  
 EMEA +31.33.45.45.600 / [info.eu@ergotron.com](mailto:info.eu@ergotron.com)  
 Custom: [custom@ergotron.com](mailto:custom@ergotron.com)

© Ergotron, Inc. 09.08.2021 Content subject to change. Patent information available at [www.ergotron.com/patents](http://www.ergotron.com/patents).  
 Ergotron devices are not intended to cure, treat, mitigate or prevent any disease.  
 WorkFit, StyleView, LearnFit, Neo-Flex, TRACE, PowerShuttle, LiFeKinnex, CareFit and eKinnex are trademarks of Ergotron.

