


ergotron®

moving you forward

StyleView 42 Medical Cart


- Easy to maneuver: Moves easily across carpet or hard surface floors with a small footprint that fits in space-constrained areas
- Enhances well-being and productivity: Caregivers can choose their preferred height for ergonomic typing and charting as they sit or stand
- Promotes patient and caregiver safety: Entire power system and cart are certified to UL/EN/IEC 60601-1 and comply with EMI limits of FCC Part 15 Class A
- Long battery life: LiFe battery option helps improve uptime with an ultrafast recharge and long runtime that can last up to 12 hours, depending on the number of connected devices
- Quick integration: Ships fully assembled



LAPTOP



LCD ARM



LCD PIVOT

VERSATILE POWERED CART SUPPORTS
MULTIPLE WORKFLOWS

STYLEVIEW™

StyleView Medical Carts bring comfort and flexibility to any workflow by adapting to meet the needs of every caregiver. Built on decades of experience, these compact carts offer top-notch ergonomics and versatility. The open architecture design, broad portfolio of accessories and quality construction makes StyleView a long-term investment that evolves with your organization.

Customize your fleet to fit your workflows with:

Three display mount options

Lithium Iron Phosphate (LiFe) or Sealed Lead Acid (SLA) battery power

Proprietary Constant Force™ (CF) manual height adjustment

STYLEVIEW™

PRODUCTS INCLUDE:

LCD Pivot or LCD Arm mounting options, antimicrobial worksurface, CPU lock, secure CPU compartment, four-port USB hub, scanner holder bracket, mouse holder, wrist re negative-tilt keyboard tray, left/right mouse tray, rear handle, storage basket, cable str four casters, eKinnex™ fleet management software.

- 1 Convenient scanner holder bracket supports smooth workflows
- 2 Large worksurface and wrist rest are treated with antimicrobial to support infection control
- 3 Negative-tilt keyboard tray, wrist rest and mouse holder support comfortable documentation; four amber LEDs provide visibility of the keyboard tray in low-light care areas
- 4 Rear storage basket keeps task-related items close for easy access
- 5 T-slot channel interface on the column allows you to integrate accessories
- 6 Four 10.2 cm lockable casters, including one with a tracking function to easily steer while moving the cart
- 7 Compact base fits into crowded healthcare spaces



LCD PIVOT CART
CONSTANT FORCE, LiFe

Durable, easy-to-clean exterior composed of aluminum, high-grade plastic and zinc-plated steel



T-SLOT WIPES HOLDER
98-450



4 DRAWERS (1 PRIMARY TRIPLE AND
1 SUPPLEMENTAL SINGLE)



9 DRAWERS

Additional accessories available at healthcare.ergotron.com. Drawers can be added or removed based on changing needs. Please note, adding drawers may impact the total height adjustment.

| | CF Laptop, LiFe | CF LCD Arm, LiFe | CF LCD Pivot, LiFe | CF Laptop, SLA | CF LCD Arm, SLA | CF LCD Pivot, SLA |
|---|--|-----------------------------------|--------------------|----------------|-----------------------------------|-------------------|
| Part # (white) | SV42-6102 | SV42-6202 | SV42-6302 | SV42-6101 | SV42-6201 | SV42-6301 |
| Typical LCD size | ≤ 17.3" | ≤ 24" | | ≤ 17.3" | ≤ 24" | |
| <small>LCD SIZE IS APPROXIMATION. DIMENSION IS DIAGONAL MEASUREMENT OF SCREEN. LCD SIZE MAY BE EXCEEDED AS LONG AS SCREEN WEIGHT DOES NOT EXCEED MAXIMUM WEIGHT CAPACITY OF THE MOUNT</small> | | | | | | |
| Weight Capacity* | ≤ 8.2 kg | ≤ 14.1 kg | | ≤ 8.2 kg | ≤ 15.9 kg | |
| Lift* | 51 cm | 68.6 cm | 63.5 cm | 51 cm | 68.6 cm | 63.5 cm |
| Tilt* | — | 65°/+60°/-5° | 25°/+20°/-5° | — | 65°/+60°/-5° | 25°/+20°/-5° |
| Pan* | — | 180° | 24° | — | 180° | 24° |
| Rotation* | — | 360° | 90° | — | 360° | 90° |
| VESA | — | MIS-D (100 x 100 mm & 75 x 75 mm) | | — | MIS-D (100 x 100 mm & 75 x 75 mm) | |
| Cart Weight | 42.3 kg | 44.1 kg | | 57 kg | 58.5 kg | |
| Shipping Dimensions | 142 x 86 x 62 cm | | | | | |
| Shipping Weight | 57.7 kg | 59.5 kg | | 72 kg | 73.5 kg | |
| Warranty | 5-year warranty for mechanical components; 3-year warranty on LiFe batteries, electric lift and switch; 2-year warranty on power modules and power coiled cords. See Ergotron product warranty terms & conditions for more detail.* | | | | | |

*See website

